

2017

Product Catalogue

Swedish premium quality

Competitive prices

Proven performance



Introducing Swedish Implant Technology

Swedish Implant Technology (SIT) was born out of a desire to provide sophisticated dental implant products at a competitive price. We are now pleased to present Osstera[®], a well proven product range of Swedish premium quality.

Our extensive industry experience has made us who we are: a specialized partner that always puts personal customer relationships and dedicated service first.

Swedish Implant Technology's management system fulfills the requirements of SS-EN ISO 13485:2012 with respect to development and sales of dental implants and accessories.

Contents


Introduction	02
Osstera [®]	04
Trilobe	10
Conical	12
Tissue level	14
Bone level	16
Contact	18



Sophisticated
technology at a
competitive price




The well proven Osstera® products are characterized by Swedish premium quality and a competitive price. The implants are designed to suit every clinical situation. The Osstera® have four different prosthetic interfaces.



Modern methods. Proven products.

Developed at the University of Gothenburg, Sweden, Osstera's patented surface consists of a nano-porous titanium oxide layer, saturated with 11% calcium. The unique combination of calcium deposits and 3-dimensional interconnecting nano-porosity makes the Osstera's surface replicate nature closely.

This surface structure has been clinically proven through 10 years of research to support successful osseointegration and has shown excellent long-term results. Numerous articles published globally also support the efficiency of Osstera's surface. To view or download research material, please visit our website www.sitdental.com



The surface has been developed at the University of Gothenburg

CAD/CAM Libraries

The Osstera CAD/CAM libraries are available for Exocad, Dental wings and 3Shape® software and can be download at our website.



Osstera® Trilobe

TRILOBE

The Osstera Trilobe with it´s tapered/parallel body is designed for high initial stability.



Trilobe Implant

With cover screw

Length mm.	8.5	10.5	11.5	13.5
NP	403508	403510	403511	403513
RP	404308	404310	404311	404313



Replica

NP	453500
RP	454300
WP	455000



Esthetic Abutment

With abutment screw and replica

NP	463533
RP	464333
WP	465033



Abutment Screw Unigrip

NP	473533
RP/WP	474333



Healing Abutments

Height mm.	3.0	3.0w	5.0	5.0w
NP	433503	433513	433505	433515
RP	434303	434313	434305	434315



Plastic Abutment

With abutment screw and replica

	Engaging	Non-Engaging
NP	463511	463510
RP	464311	464310
WP	465011	465010



Titanium Base

With abutment screw and replica

	Non-Eng. 0.4 mm	Eng. 1.0 mm	Non-Eng. 1.0 mm
NP	463538-00	463538-11	463538-10
RP	464338-00	464338-11	464338-10



Multi Unit Abutment

Height mm	1.0	2.0	3.0
NP	463539-1	463539-2	463539-3
RP	464339-1	464339-2	464339-3

Height mm	2.0 17°	3.0 17°	4.0 17°	4.0 30°
NP	463539-217	463539-317		
RP	464339-217	464339-317	464339-417	464339-430



Impression Coping

NP	443500
RP	444300



CoCr Abutment

With abutment screw and replica

	Engaging	Non-Engaging
NP	463521	463520
RP	464321	464320
WP	465021	465020



Scanning Abutment

NP	493511
RP	494311




Healing Cap	430010
Impression Coping	440030
Replica	450000
Plastic Cylinder*	460010
CoCr Base*	460020
Titanium Base*	460030
Abutment Screw	470033
Scanning Abutment	490011

*With abutment screw and replica

Osstera® Conical

CONICAL

The Osstera Conical combines the Trilobe implant body with a conical prosthetic connection.



Conical implant

With cover screw

Length mm	8.5	10.5	11.5	13.5
NP	503508	503510	503511	503513
RP	504308	504310	504311	504313



Replica

NP	553500
RP	554300



Esthetic Abutment

With abutment screw and replica

NP	563533
RP	564333

Abutment Screw Unigrip

NP	573533
RP	574333



Healing Abutments

Height mm	3.0	3.0w	5.0	5.0w
NP	533503	533513	533505	533515
RP	534303	534313	534305	534315

Plastic Abutment

With abutment screw and replica

	Engaging	Non-Engaging
NP	563511	563510
RP	564311	564310



Titanium Base

With abutment screw and replica

	Non-Eng. 0.4 mm	Eng. 1.0 mm	Non-Eng. 1.0 mm
NP	563538-00	563538-11	563538-10
RP	564338-00	564338-11	564338-10

Impression Coping

NP	543500
RP	544300

CoCr Abutment

With abutment screw and replica

	Engaging	Non-Engaging
NP	563521	563520
RP	564321	564320



Scanning Abutment

NP	593511
RP	594311



Osstera® Tissue level

TISSUE LEVEL

Osstera Tissue Level implants are designed for one-stage surgical procedures. It saves time and increases efficiency.



Tissue level Implants RP

With cover screw

Length mm	8.0	10.0	12.0	14.0
Ø3.3	303308	303310	303312	303314
Ø4.1	304108	304110	304112	304114
Ø4.8	304808	304810	304812	304814



Replica

RP 354800



Synocta Abutment

RP 364833



Scanning Abutment

RP 394811



Healing Abutments

Height mm	2.0	3.0
RP	334802	334803



Plastic Abutment

With abutment screw and replica

	Engaging	Non-Engaging
RP	364811	364810



Synocta Abutment Angulated

	15°A	15°B	20°A	20°B
RP	364834-A	364834-B	364835-A	364835-B

A = angle to the edge of the Octagon
B = angle to the flat surface of the Octagon

Abutment Screw Torx

RP 374832



Impression Coping

RP 344800



CoCr Abutment

With abutment screw and replica

	Engaging	Non-Engaging
RP	364821	364820



Titan Base

With abutment screw and replica


	Non-Eng. 0.4 mm	Eng. 1.0 mm
RP	364838-00	364838-11



Osstera® Bone level

BONE LEVEL

The Osstera Bone Level implant stands for excellent clinical performance and patient satisfaction.



Bone level Implants

With cover screw

Length mm	8.0	10.0	12.0	14.0
Ø3.3	603308	603310	603312	603314
Ø4.1	604108	604110	604112	604114
Ø4.8	604808	604810	604812	604814



Replica

NP	653300
RP	654100



Titanium Base

With abutment screw and replica

	Non-Eng. 1.0 mm	Eng. 1.0 mm
NP	663338-00	663338-11
RP	664138-00	664138-11

Healing Abutments

Height mm	2.0	4.0	6.0
NP	633302	633304	633306
RP	634102	634104	434106

Plastic Abutment

With abutment screw and replica

	Engaging	Non-Engaging
NP	663311	663310
RP	664111	664110



Scanning Abutment

NP	693311
RP	694111

Impression Coping

NP	643300
RP	644100

CoCr Abutment

With abutment screw and replica

	Engaging	Non-Engaging
NP	663321	663320
RP	664121	664120



Abutment Screw Torx

NP/RP	673332
-------	--------



Contact us

To get more information about our products and updates,
or to make a purchase, please visit www.sitdental.com.
We are also available at +46 31 29 66 00.



